



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Classical Mechanical Pipette

Suitable for pharmaceutical and chemical companies, school laboratories etc.

FEATURES:

- Light weight and ergonomic design
- Digital display, easy to read data
- Volume range: 0.5 μ l to 1ml
- Easy to calibrate and maintain
- The lower part can be autoclaved without affecting performance
- Manufactured from innovative materials



0.5-10 μ l

2-20 μ l

10-100 μ l

20-200 μ l

100-1000 μ l

MODEL SPECIFICATIONS:

Volume Range	Increment	Test Volume (μl)	Capacity Tolerance (%)	Test Repeatability (%)
0.5-10 μl	0.1 μl	1	± 2	≤ 3
		5	± 1.5	≤ 2
		10	± 1	≤ 2
2-20 μl	0.1 μl	2	± 2	≤ 3
		10	± 1.5	≤ 2
		20	± 1	≤ 2
10-100 μl	1 μl	10	± 2	≤ 2
		50	± 1.5	≤ 1
		100	± 1	≤ 1
20-200 μl	1 μl	20	± 2	≤ 2
		100	± 1.5	≤ 1
		200	± 1	≤ 1
100-1000 μl	5 μl	100	± 2	≤ 1
		500	± 1	≤ 0.5
		1000	± 1	≤ 0.5