



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein
Pretoria, Gauteng, 0081


ANALYTICAL EXCELLENCE...

www.microanalytica.co.za

Forced Air Drying Oven MA-VF-Series

FEATURES:

- External material: cold rolled steel with anti-bacteria powder coating
- Stainless steel inner chamber, with adjustable shelves
- PID microprocessor intelligent temperature controller
- Double-layered glass observation window, adjustable airtight buckle lock (except MA-V35F)
- Over-temperature protection
- Optional: Independent temperature protector





MA-V35F



MA-V70F/V136F/V225F/V429F/V640F/V960F



MA-V30F/V45F/V65F/V125F



MA-V230F



MA-V625F

MODEL SPECIFICATIONS:

Model	MA-V30F	MA-V35F	MA-V45F	MA-V65F	MA-V70F	MA-V125F
Capacity	30L	35L	45L	65L	70L	125L
Temp. Range	RT+10~300°C	RT+10~250/300°C	RT+10~300°C	RT+10~300°C	RT+10~250/300°C	RT+10~300°C
Temp. Precision	0.1°C					
Temp. Fluctuation	±0.5°C					
Ambient Temp.	5~40°C					
Timing Range	1~9999min					
Shelves No.	2 pcs					
Display	LED	LCD	LED		LCD	LED
Consumption	800W	550W	1200W	1600W	1250W	2300W
Power Supply	AC220V±10%, 50/60Hz (standard); AC110V±10%, 60Hz (optional)					
Internal Size (W*D*H)	310*310*310	340*325*290	350*350*350	400*360*450	400*375*470	500*450*550
External Size (W*D*H)	460*510*695	490*520*690	500*550*735	550*550*840	540*550*870	636*680*915
Package Size (W*D*H)	550*580*770	620*630*830	620*585*805	640*640*905	670*660*1040	725*735*1005
Gross Weight (kg)	34	55	39	47	75	58

Model	MA-V136F	MA-V225F	MA-V230F	MA-V429F	MA-V625F	MA-V640F	MA-V960F
Capacity	140L	225L	230L	420L	625L	600L	900L
Temp. Range	RT +10~ 250/300°C	RT +10~ 250/300°C	RT +10~300°C	RT +10~ 250/300°C	RT +10~300°C	RT +10~250/300°C	RT +10~ 250/300°C
Temp. Precision	0.1°C						
Temp. Fluctuation	±0.5°C						
Ambient Temp.	5~40°C						
Timing Range	1~9999min						
Shelves No.	2 pcs			3 pcs		4 pcs	
Display	LCD		LED	LCD	LED	LCD	
Consumption	1500W	2100W	3000W	4000W	5000W	4500W	6000W
Power Supply	AC220V±10%, 50/60Hz (standard); AC110V±10%, 60Hz (optional), 380V, 50Hz (MA-V225F - 300°C)			Standard: 380V, 50/60Hz Optional: 220V, 50/60Hz	380V, 50/60Hz		
Internal Size (W*D*H) mm	450*557*520	500*615*720	600*500*750	600*550*1300	660*760*1250	800*600*1300	1000*600*1600
External Size (W*D*H) mm	600*740*920	650*800*1130	730*670*1220	750*730*1700	895*925*1830	950*780*1720	1140*800*2150
Package Size (W*D*H) mm	730*860*1080	780*920*1270	770*840*1420	870*860*1830	995*980*1980	1080*910*1880	1260*910*2200
Gross Weight (kg)	85	110	103.5	145	277	170	230