



+27 60 999 9140

info@microanalytica.co.za

489 Jacqueline Drive, Garsfontein  
Pretoria, Gauteng, 0081

ANALYTICAL EXCELLENCE...

[www.microanalytica.co.za](http://www.microanalytica.co.za)

### **Multifunctional Incubator MA-200 / MA-300**

Widely used in the storage of medicine, reagents, vaccines, dialysate heating, heat preservation and preservation of medical equipment in the operating room, normal saline 37°C – 42 °C heating and washing fluid 35°C - 37°C heating etc.

#### **FEATURES:**

- Unique air duct design
- Non-slip separator design prevents the sample from falling off when the shelf is pulled
- 7-inch colour LCD touch screen can monitor the temperature curve in real time, menu-type operation interface is easy to understand
- 3-level operation interface authority (administrator, experimenter, operator)
- Equipped with an interconnected PC interface for remote operation monitoring and optional mobile phone APP control program




MA-200




MA-300



LCD display



UV & LED lamp



USB LAN port

**MODEL SPECIFICATIONS:**

<b>Model</b>	<b>MA-200</b>	<b>MA-300</b>
Capacity	200L	300L
Temp. Range	0~100°C	
Control Type	Microcomputer Controller	
Temp. Fluctuation	≤±1°C (at 37°C)	
Temp. Uniformity	≤±2°C (at 37°C)	
Time Setting	1-5999min	
Alarm	Audible and visual alarm; over-temperature alarm	
Shelves	3 pcs/adjustable	
Material	Exterior: cold-rolled steel with painting; interior: stainless steel	
Door	Double tempered glass doors with heating function	
Sterilization Type	UV sterilisation lamp	
Ambient Temp & Humidity	15~30°C, suggest 25±2°C; relative humidity ≤80%	
Caster	Universal caster	
Consumption	1200W	1450W
Power Supply	AC220V±10%, 50/60Hz (standard); AC110V±10%, 60Hz (optional)	
Internal Size (W*D*H) mm	518*500*832	500*490*1295
External Size (W*D*H) mm	598*797*1475	761*606*1900
Package Size (W*D*H) mm	686*885*1618	801*646*2020
Gross Weight (kg)	130	156